


Cyber Ho 60

HOLMIUM LASER





Cyber Ho 60

THE REVOLUTION IN LITHOTRIPSY

Cyber Ho Holmium laser (2.1 μm) meets the increasing demand of efficacy, flexibility and cost-effectiveness with a unique multi-application laser platform able to perform both **Lithotripsy** and **HoLEP**.


Cyber Ho brings outstanding innovation by offering the exclusive **Vapor Tunnel™**, **Virtual Basket™** and **MasterPULSE** technology for advanced retropulsion control, along with impressive energy and frequency (up to **60 Hz**) settings.

General Overview

- ✓ Effective Lithotripsy
- ✓ High Frequency Emission (up to 60 Hz)
- ✓ Minimized Retropulsion
- ✓ BPH Treatment
- ✓ Reduced Depth of Penetration (0.3 - 0.4 mm)
- ✓ Soft Tissue Surgery
- ✓ High Versatility
- ✓ Quick ROI




Double Footswitch



The double footswitch enables **immediate** switch from one emission mode to another, with **complete customization** of pedal-mode association. No bothersome interruptions are needed for settings readjustment.

Fragmentation



**SHORT PULSE
(HIGH ENERGY)**



HIGH PULSE ENERGY

5 J for strong pulse effect



TREAT ALSO THE HARDEST STONES

Greater pulse energy allows to break harder stones



COLLECTION BASKET NEEDED

Retrieve stone pieces upon fragmentation

Dusting Effect



LIMITED RETROPULSION

Easy ablation with no need to chase the stone



NO NEED FOR BASKET

The obtained fine dust obviates the retrieval phase



LONG PULSE WIDTH

Up to 1100 μ s, for smooth *Long Pulse Dusting*



EXTREME FREQUENCY

Up to 60 Hz, for enhanced speed in *Dusting* action



**HIGH
FREQUENCY**

Vapor Tunnel™

Consisting in a Single **Specific Long Pulse**, this emission mode allows limited retropulsion and fine stone ablation. The Vapor Tunnel™ is designed in order to use the minimum peak power in accordance with selected output settings.




Bubble Dynamics of Vapor Tunnel™



Virtual Basket™

Composed by a **Double Pulse Emission***, combines a low retropulsion with a fragment suction effect.

*(time duration separating the two pulses is chosen so that the second pulse is emitted from the distal tip of the fiber when the bubble size, and the corresponding amount of displaced fluid, is at a maximum)



Bubble Dynamics of Virtual Basket™



Advantages of Virtual Basket™ & Vapor Tunnel™



NO EXTRA COSTS

These modes do not need dedicated and more expensive fibers, bringing the mentioned advantages without extra expenses.



MAGNETIC EFFECT

These modes allow stone ablation while holding the target in place, without inducing stone retropulsion.



EASIER TREATMENT


With a more stable target, lithotripsy treatment can proceed easily with fewer hassles.



TIME SAVING

Less stone retropulsion prevents the time-consuming fiber repositioning, whereas enhanced energy transmission increases the ablation rate.


This long bubble touching the target represents a direct connection between fiber tip and stone, granting enhanced energy delivery



As the pulse ends, the bubble collapses.
The stone is dragged backwards together with the collapsing bubble (like a virtual basket)



MasterPULSE




General Benefits

Reducing retropulsion and modifying the tissue cutting get easier: instead of trying multiple different settings, start with your preferred ones and then adjust the MasterPULSE to tune the effect of laser emission based on your visual feedback. Regulation of pulse width has never been so easy!



EASE OF TREATMENT

Experience a more intuitive and different way to adjust laser settings, simply based on your visual feedback



Example of Lithotripsy mode



REDUCED STONE INSTABILITY

Lower the stone instability step by step, by progressively increasing pulse width



CUTTING DOWN TREATMENT TIME

Obtain the desired effect quickly, without losing time with the standard adjustment of energy and frequency parameters



GREATER FLEXIBILITY

7 levels of pulse width offer a greater flexibility with respect to the traditional 3 levels offered by the other holmium devices

BPH

The large amount of literature demonstrates the efficacy and safety advantages of **HoLEP** (*Holmium Laser Enucleation of the Prostate*) with respect to traditional treatments available for BPH (*Benign Prostatic Hyperplasia*). Recent studies and trials have validated the excellent outcomes achieved by this technique, with its success being reproduced in a diverse array of patients.

Clinical evidences proved a 50 - 60 W Holmium output to work perfectly for an effective and gentle enucleation. HoLEP with Cyber Ho provides excellent results in terms of **long-term effectiveness, safety, hemostasis and fast recovery**.



Starting 5 o'clock incision in HoLEP

- ✓ Precise cutting
- ✓ User preset saving
- ✓ Guided selection of settings
- ✓ Double footswitch



Quick switch from one emission mode to another (eg. from cutting to coagulation emission)



Delicate Areas Mode

Dedicated Coagulation Setting

User-friendly Software



GUIDED SELECTION

- BPH
- Lithotripsy
- Soft Tissues

SAVE AND LOAD OF SETTINGS



Fibers

Cyber Ho device can be operated with a large range of fibers, depending on the application, flexibility and settings required.



STANDARD FIBERS

For general use in stone and soft tissue treatments



BALL TIP FIBERS

Strongly simplify the insertion in already bent scopes



PERFORMANCE FIBERS

Designed to deliver great power even with small fibers



GASTRO FIBERS

Specifically designed for the fragmentation of gallstones



FIBER RECOGNITION

Cyber Ho automatically adjusts emission settings based on the connected fiber diameter



AVAILABLE DIAMETERS

200, 272, 365, 550, 800 and 1000 μm



REUSABILITY

All fibers are available both as disposable and reusable (except ball tip model)



CLEANING

Reusable fibers can be sterilized by Sterrad[®] and steam sterilization




STERILIZATION TRAY

A dedicated tray for sterilization of fibers and tools

Applications

Cyber Ho 60 can be used to perform incision, excision, resection, ablation, vaporization, coagulation and hemostasis of soft tissue and in lithotripsy of stones in various medical specialties, for example:



Technical Specifications

Wavelength	2,1 μm
Average power	Up to 60 W
Repetition rate	Up to 60 Hz
Energy per pulse	Up to 5 J
Pulse duration	50 ÷ 1100 μs
Beam delivery	Wide range of flexible silica fibers
Aiming beam	532 nm (adjustable <5 mW) - Class 3R
Fiber recognition	RFID System
Activation	Double footswitch
Electrical requirements	200-230 Vac; 50/60 Hz; 5 kVA, single phase
Cooling	Internal chiller
Operating temperature	10°C ÷ 30°C
Laser class	4
Dimensions and weight	52 cm (W) x 120 cm (D) x 123 cm (H) (monitor closed) - 200 kg

*This brochure is not intended for the U.S. market.
Certain Intended Uses/Configurations/Models/Accessories are not cleared for U.S.*

VISIBLE AND INVISIBLE LASER RADIATION

Avoid eye skin exposure to direct or scattered radiation

Laser product: Class 4

Aiming beam: Class 3R



CAUTION - Laser radiation



Note: National local authorities may put restrictions to the parameters indicated in the table in the previous page, or may limit or remove certain intended uses. Specifications are subject to change without notice.

Quanta System products are manufactured according to the International standards and have been cleared by the most important International notified bodies.

The Company is UNI EN ISO 9001:2015 and EN ISO 13485:2016 certified. Quanta System S.p.A. was founded in 1985 and belongs to the El. En. Group (a public company listed in the Star segment of the Italian Stock Exchange) since January 2004.

The company, divided into three business units (medical, scientific and industrial) is specialized in manufacturing of laser and opto-electronic devices.

Quanta System
LASER IN OUR DNA

