

TryTable

Multipurpose examination couch – solution for urology and gynaecology, USG scanning, X-ray images




↑
"Puzzle" – easy to release / fix

- Wide range of settings
- Wide range of accessories
- 3 motor examination couch
- Universal eurolist to attach the accessories
- Easy to clean/disinfect
- Four locking wheels
- Levers easy to use
- Accessories right/left side



0°–45°



Headrest easily adjustable



3 motor examination couch / Memory for 3 positions



0°–85°

Back angle



650–950 mm

Height



0°–15°

Seat angle



Accessories

Wide range of supports for urology and gynaecology examinations



↑
Hand rests



↑
Leg supports



↑
Leg supports



↑
Leg rests

Universal holders for easy fixing / releasing



↑
Lever I.



↑
Lever II.



↑
Stainless bowl



↑
Foot pedal



TECHNICAL SPECIFICATION

- Working area width 590 mm
- Operational temperature +15°C – +35°C
- Working area length 1.400 – 1.900 mm
- Power supply 230V, 50Hz
- Total permitted loading 150 kg
- Maximum power input 500 VA
- Weight 120 kg