



Uromic Jive

System which always adapts
to your daily needs.

- Compact 5 channel system for complex functional diagnostics.
- Height adjustable stand.
- Device chassis equipped with braked wheels allows for easy movement of the device, facilitated by convenient handrails.
- Small footprint minimises use of space in clinic and storage.
- Infection control friendly keyboard.
- Unique wireless uroflowmeter connection at 866/915 MHz ensures the stability and no signal dropout.

Various options upon customer's choice:


- wired/wireless keyboard with touchpad connection,
- wide shelf for standard or medical keyboard and mouse,
- widescreen, standard or touch monitor.



Version with a shelf with
medical waterproof white
mouse and keyboard.



Version with a shelf with
standard mouse and
keyboard.



Easy to release/fix the
height with one screw.