



Endoscopy LED Light source

U-Lite 200LED

Special Feature:

- Bright White day light economical LED Light Source
- Easy to use with standard fibre optic connections.
- Inbuilt Temperature sensor to avoid overheating
- Light weight & long-life LED Lamp
- Less fan noise with high cooling effect.
- Digital Light intensity controller with display
- Extremely high Optical Output from a 14mm² circular source.

Connector for Fiber Optic Cable



Technical Data:

- Lamp : Monolith Large Chip LED Technology
- Colour Temperature : Up to 6500K
- Lifetime : Up to 50,000 Hours
- Colour : Daylight White
- Current Rating : Min. Drive current 0.2Amp to Maximum current
28Amp
- Luminous Flux : 2400 – 2600 Lumen
- Maximum Junction Temperature : 150° C
- Storage Temp. Range : -40/ +100° C

Electrical System

- Power connections: 220-240VAC/ 50-60Hz/ 1.0A-5.0Amp
- Power consumption: 150W

UNIQUE
Medical Devices

324, Vipul Business Park, Sector -48,
Gurgaon – 122001, Haryana. INDIA
Email- sales@uniquemedicaldevices.org & info.umd@gmail.com
[Mob. -0124-4080411/ 8860086235/ 8860086234](tel:-0124-4080411)
[Web- www.uniquemedicaldevices.org](http://www.uniquemedicaldevices.org)