



Endoscopic Camera



1-Chip CMOS HD

UMD-UCAM-21

FUNCTION

HIGH RESOLUTION

Output video up to 1920*1080P60fps

RICH OUTPUT INTERFACE

Includes HDMI(2pcs),DVI,SDI,CVBS,USB

BUILT-IN HDMI OUTPUT

The processing board built-in a extra HDMI interface to realize the camera + monitor integrated machine

FRIENDLY SPECIAL FUNCTION

User defined:
Color, Surgical scenes, Shortcut Buttons...

PLAYBACK VIDEO

Supports recording 1080p60fps video ,
Supports playback recorded video/photos.

SUPPORTS RS232 TOUCH MENU

You can make the touch menu LCD on the front panel of the machine, instead of buttons

4 BUTTONS ON CAMERA HEAD

The camera module has up to 4 buttons,
Each button can be user defined, Like

AWB,FREEZE,BRIGHTNESS,ZOOM,RECORD,PHOTO





Technical Specificaation

- Camera Resolution : 1920 X 1080 P60FPS
- Sensor (Chip) : 1/1.8" CMOS Technology
- Resolution : " >1000 TV Line
- Noise Level : 55 dB
- Sensitivit y : 2000 1XF11/ 3200K
- Surgical Scene : Laparoscopy ENT, Hyteroscopy, Urology, Fiberscope
- Recording Mode : MP4
- Record time : 30 min./file
- Photo : JPEG
- Output Mode : SDI, CVBS, DVI, HDMI


Operating Conditions :

- Temperature : 10- 40° C
- Output Mode : SDI, CVBS, DVI, HDMI

Input Electrical Rating

- Frequency : 50Hz
- Voltge / Current : 220VAC/ 5Amp





UNIQUE
Medical Devices

324, Vipul Business Park, Sector -48,
Gurgaon – 122001, Haryana. INDIA
Email- sales@uniquemedicaldevices.org & info.umd@gmail.com
[Mob. -0124-4080411/ 8860086235/ 8860086234](tel:-0124-4080411)
[Web- www.uniquemedicaldevices.org](http://www.uniquemedicaldevices.org)