

# Deepra-Flex - Slim

CMOS - Flexible VIDEO Ureteroscope



**UNIQUE**  
Medical Devices

# Deepra-Flex Slim (8.5 Fr. Video URS)




Video System+Light Source+Monitor



Rotation Function, Left:120° Right:120°



Outer diameter 3.0mm  
with 1.2mm instrument channel



## Insertion Rotatable Part

- Smooth Insertion
- 120 Degree left/ right Movement
- Improved Therapeutic Capability
- Easier Access

## Maximum Deflection

- 275 Degree Both Side
- Logical & Anti-Logical Movements available
- Improved shaft stiffness for better steering to allow to reach lower clayx easily

Optical system	Field of View	120°
	Direction of View	Forwarding
	Depth of Field	2-50mm
Insertion Tube	Distal End Outer Diameter	8.5 Fr.
	Inner Channel Diameter	3.6Fr(1.2mm)
	Flexible Tube Outer Diameter	8.5 Fr.
	Working Length	670mm
	Angulation Range	UP 275° /DOWN 275°
Monitor	Power	DC12V
	Battery Working Time	4 hours
	Screen Size	8 inch
	Resolution	680 x 480
	Recording Function	Yes
	Storage	8GB
Optional	Trolley	1.25 meter

